
INSTI[®] HIV SELF TEST INSTRUCTIONS

Intended Use:
for untrained lay users as a self test to detect HIV-1 and HIV-2 antibodies using a drop of human fingerstick blood


Questions?
+1-604-204-6784
www.INSTI-HIVSelfTest.com


Do not use if the test device pouch is broken.

Do not use if you:
• have a bleeding disorder
• are on ART

PREPARATION

STEP 1: COLLECT BLOOD

STEP 2: TEST

STEP 3: READ RESULT


READ RESULT RIGHT AWAY AND **WITHIN 5 MINUTES.**

NEGATIVE

3 TEST AGAIN IN 3 MONTHS

INSTI[®] has a specificity of 99.5%. This means a negative result will be correct 995 out of every 1000 tests.

POSITIVE

It does not matter how light or dark either dot is.

GO TO CLINIC FOR CONFIRMATORY TESTING

INSTI[®] has a sensitivity of 99.8%. This means a positive result will be correct 998 out of every 1000 tests.

INVALID

GO TO CLINIC FOR CONFIRMATORY TESTING

DISPOSAL

1. Put all items back into the pouch. 2. Throw away pouch in waste bin.

Blood can transmit infectious diseases. Clean up spills.

Dispose in accordance to local regulations.

Manufacturer Caution Use by date LOT Lot number
REF Catalogue number IVD For in vitro diagnostic use only Do not reuse X_n Contains sodium azide Harmful if swallowed
Consult the Self Test instructions Store at 2 - 30°C or refrigerated (2 - 8°C) if required Sterilization using irradiation

Manufactured by bioLytical Laboratories Inc.
1108 - 13351 Commerce Parkway,
Richmond, BC V6V 2X7 Canada
Phone: +1 604 204 6784
Fax: +1 604 244 8399
www.bioLytical.com
©Copyright 2017 All rights reserved.

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Imetengenezwa na bioLytical Laboratories Inc. 1108 - 13351 Commerce Parkway, Richmond, BC V6V 2X7 Kanada Simu: +1 604 204 6784 Faksi: +1 604 244 8399 www.biolytical.com ©Copyright 2017 Haki zote zimehifadhiwa. 51-1192(D) 01-May-2017